## ENGLISH

These instructions were developed using the guidance given in the following standards:

- ISO 17664: Processing of health care products – Information to be provided by the medical manufacturer for the processing of medical devices.
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.
- ANSI/AAMI ST179: Comprehensive Guide to Steam Sterilization and Sterility Assurance in health care Facilities.
- ISO 16061: Instrumentation for use in association with non-active surgical implants – general requirements.
- Guidance for Industry and Food and Drug Administration Staff: “Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling”. June 9, 2017.

#### A. INSTRUMENT SET INFORMATION AND DESCRIPTION

Intended use:

The WasherCap™ Surgical Instrument Set is a reusable surgical instrumentation intended to be used in the context of a planned Anterior Cruciate Ligament (ACL), carried out with the WasherCap™ Fixation System with the aim of ensuring that the surgical procedure and the use of the WasherCap™ are performed correctly. The WasherCap™ Surgical Instrument Set is composed of reusable surgical instruments, made of surgical stainless steel according to ASTM F-899 and ASTM F138. Each Set is composed of:

- 1 Sterilization tray.
- 4 Countersink Positioners for the complete range of the WasherCap™ device (7mm, 8mm, 9mm and 10mm).
- 1 Positioner guide.
- 2 Kirschner wires.
- 4 Countersinks (consisting of 2 pieces: 1 Main body and 1 Stop element) for the complete range of the WasherCap™ device (7mm, 8mm, 9mm and 10mm).
- 1 Handle with AO connection system.
- 2 Inserters for the complete range of the WasherCap™ device (7mm & 8mm and 9mm & 10mm).
- 2 Inserter Coupling Screws for the complete range of the WasherCap™ device (7mm & 8mm and 9mm & 10mm).
- 1 Punch.
- 2 Aaws for the complete range of the WasherCap™ device (7mm & 8mm and 9mm & 10mm).
- 1 WasherCap™ Screwdriver.
- 1 Countersink tool.

Abanza provides detailed surgical techniques in print, video and electronic formats. Abanza also provides demonstrations and information on detailed surgical techniques (www.abanzamed.com).

#### B. LIMITATIONS ON REPROCESSING

Repeated reprocessing has minimal effect on these reusable instruments when performed according to ABANZA's instructions. End of useful life is normally determined by wear and damage due to use. In addition, it is recommended that, prior to re-using the instruments (after cleaning process), a visual inspection is carried out to verify that no signs of visible dirt, corrosion and/or other type of damage are present. Please contact your local distributor or Abanza representative for more information on the reuse test and if you detect wear and/or damage on the instruments for replacement.

#### C. REPROCESSING VALIDATION

When properly performed cleaning, disinfection and sterilization, do not compromise the use and mechanical performance of these instruments. The instruments are used with or on patients who may be harboring either recognized and unrecognized infections. To prevent the spread of infection, all reusable instruments must be thoroughly cleaned, disinfected, and sterilized after use in each patient, following the specifications of this Instructions for Use.

The recommended cleaning, disinfection and sterilization methods in this IFU have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the half cycle "overkill" approach was used for sterilization validation and demonstrates a sterility assurance level (SAL) of 10<sup>-4</sup>. Abanza has validated the reprocessing instructions described in this document and, in more detail, in the document "Instructions for care, cleaning, maintenance and sterilization" (CI-001) as capable and effective for the reuse of the instruments which compound the set. However, it remains the responsibility of the responsible user to ensure that the reprocessing of the instruments is carried out correctly and effectively, in compliance with legislation and ordinances and using equipment and material approved by the FDA, CE marked or authorized by the local government body responsible for them.

#### D. CLEANING AND DISINFECTION

- CLEANING AT THE POINT OF USE AND HANDLING IN THE REPROCESSING AREA**
  - It is recommended to reprocess the instruments as soon as possible after use.
  - After use, remove excess dirt (liquid, tissue, bone fragments, organic debris, etc.) from the instruments with disposable cloth.
  - Do not allow blood, liquids or other organic debris to dry on the instruments. It is recommended that water be allowed to run through the instruments.
  - Take the instruments to the reprocessing area as soon as possible. Avoid bumping the instruments against each other during transport. It is recommended that the instruments are transported in the instruments tray, fastened to the fastening elements to avoid impacts.
  - The responsible user has the option of manual cleaning and disinfection (section D-II) or machine cleaning (automatic) and thermal disinfection (section D-III).

#### II. MANUAL CLEANING AND DESINFECTION

Manual cleaning and disinfection instructions may be used as an alternative cleaning method to the machine (automatic) cleaning and thermal disinfection.

Instruments that require disassembly should be disassembled prior to manual cleaning and disinfection. Refer to the "Instructions for care, cleaning, maintenance and sterilization" (CI-001) supplied with the instrument set to check instruments that require disassembly.

- Manual Cleaning**
  - If it has not been done previously at the point of use of the instruments, remove excess dirt from the instruments with a disposable cloth.
  - Immerse the instrument in an enzymatic or alkaline cleaning detergent solution developed specifically for this function and with proven effectiveness (FDA approved, or CE marked). Cleaning solutions should always be mixed to the manufacturer’s specification for concentration and cleaning should be conducted at ambient temperature unless otherwise stated in the cleaning solution manufacturer’s instructions.
  - Scrub the instrument with a soft bristled brush, paying special attention to areas where debris might accumulate (threads, lumens, etc). Always avoid any harsh materials that can scratch or mar the surface of the instrument.
  - Rinse the instrument thoroughly under running water until all visible traces of detergent have been removed. Actively and thoroughly wash lumens, orifices and other hard-to-reach areas.
  - Immerse the instrument in an ultrasonic bath in a neutral cleaning solution suitable for ultrasonic baths developed specifically for this function and with proven effectiveness (FDA approved, or CE marked). Prepare the solution according to the manufacturer’s instructions. Ultrasonic cleaning for a minimum of 10 minutes at 40-50 kHz is recommended, provided that the manufacturer’s instructions advise this.
  - Rinse the instrument thoroughly with purified, ultrapure or sterile water until all visible traces of detergent have been removed. Actively and thoroughly wash lumens, orifices and other hard-to-reach areas.
  - Dry the instrument with a clean cloth made of a material that does not shed any fibers. Pressurized air can be used to assist in drying.
- Disinfection (Optional)**
  - Immerse instrument in disinfection solutions. Disinfection solutions should always be mixed to the manufacturer’s specification for concentration. It is recommended to follow the manufacturer’s instructions for the correct exposure time values.
  - Remove the instruments from the disinfectant solution and rinse them with purified, ultrapurified or sterile water until all visible traces of disinfectant are removed.
  - Dry the instrument with a clean cloth made of a material that does not shed any fibers. Pressurized air can be used to assist in drying.
- Inspection after cleaning and disinfection**
  - Visually inspect each instrument and check for signs of visible dirt and/or corrosion. Pay attention to those areas where dirt may become embedded (cut areas or nooks and crannies).
  - If you detect any traces of dirt, repeat the manual cleaning and disinfection process focusing on the areas that are dirty after the process.

#### III. MACHINE (AUTOMATIC) CEANING, AND THERMAL DESINFECTION

Instruments that require disassembly must be disassembled before proceeding with the automatic cleaning and disinfection. Refer to the "Instructions for care, cleaning, maintenance and sterilization" (CI-001) supplied with the instrument set to check instruments that require disassembly.

For automatic cleaning and disinfection of instruments, washer-disinfectors should be used whose effectiveness is approved (FDA approved or CE marked according to ISO 15883). Follow the operating instructions of the washer-disinfecter manufacturer. Use cleaning agents recommended by the manufacturer or intended for use in washer-disinfectors (FDA-approved or CE marked).

- Load the instruments in the washer so that all surfaces of the devices are accessible to cleaning and any that may retain liquid can drain (hinges should be open and cannulations/holes positioned to drain).
- Run the standard automatic washing cycle for surgical instruments. The parameters of the cycle validated by Abanza are listed below:

PARAMETERS OF THE AUTOMATIC WASHING CYCLE			
Phase	Temperature	Duration	Concentration detergent
Prewashing I	10 <span> </span> °C (50 <span> </span> °F)	10 minutes	NA
Washing	55 <span> </span> °C (131 <span> </span> °F)	5 minutes	0.3 <span> </span> % detergent
Neutralization	10 <span> </span> °C (50 <span> </span> °F)	2 minutes	NA
Rinsing II	10 <span> </span> °C (50 <span> </span> °F)	1 minutes	NA

- If the surgical instrument washing cycle has a thermal disinfection part, it could be performed within the parameters set by the washer-disinfecter manufacturer, but it is optional.

- After cleaning, remove the instruments from the washer-disinfecter and visually inspect each instrument and check for signs of visible dirt and/or corrosion. Pay attention to those areas where dirt may become embedded (cut areas or nooks and crannies).
- If you detect any traces of dirt, repeat the manual cleaning and disinfection process.
- If there is excess moisture after the process, dry the instruments with a clean cloth of a non-fibrous material until they are completely dry.

#### E. MAINTENANCE AND INSPECTION

For proper maintenance, instruments should be inspected after cleaning and disinfection for damage before sterilization.

As the instruments can be affected by normal use, check that they are not damaged to the point of malfunctioning or forming projections that can damage tissue or puncture surgical gloves.

If damage is detected, do not use the device and contact your local distributor or Abanza representative for replacement.

#### F. STERILIZATION TRAY

Cleaned, disinfected and inspected instruments should be placed in the clean sterilization tray of the WasherCap™ Surgical Instrument Set. It is recommended that the instruments be placed in the fixation elements presented in the sterilization tray, specially configured for the placement of each instrument that compose the set, so that they are placed in such a way as to avoid knocks and collisions between the instruments and are arranged for optimum sterilization.

The sterilization tray should be placed in approved sterilization container at the user’s convenience.

#### G. STERILIZATION

Abanza's WasherCap™ Instrument Set is supplied non-sterile (see package labeling for more information). Before using any instrument that is not supplied in a terminal sterilization configuration, clean and disinfect the instrument and then sterilize it (see cleaning specifications described in section D of this document)

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions. Refer to "Instructions for care, cleaning, maintenance and sterilization" (CI-001) supplied with the instrument set for reprocessing.

Sterilize the assembled instruments (those that have required disassembly for the previous steps) placed on the sterilization tray and in a suitable sterilization container. In annex II of the "Instructions for care, cleaning, maintenance and sterilization" (CI-001) those cases are indicated in which it is recommended to dismantle the surgical instruments before cleaning or sterilization.

It is recommended to sterilize the instruments by moist heat method in an autoclave (sterilizer) with a pre-vacuum cycle that meets the requirements of EN 285, EN 13060, ISO 17665 and AINSI/AAMI ST179, and to perform validation and maintenance according to the same standards. Always follow the autoclave manufacturer’s instructions and check that the maximum load indicated by the manufacturer is not exceeded.

**Note:** The autoclave must be properly installed, maintained, and calibrated. The autoclave manufacturer’s operating instructions and recommended guidelines for maximum sterilization load should be followed. Your institution is responsible for the sterility assurance validation.

Abanza has validated and recommends the sterilization of the instruments according to the following parameters:

MOIST HEAT STERILIZATION PARAMETERS:

STERILIZATION PARAMETERS: FOR THE USA ONLY				
Method	Sterilizer cycle	Temperature of exposure	Exposure time	Minimum drying time
Moist heat sterilization	Pre-vacuum	132 <span> </span> °C (270 <span> </span> °F)	4 minutes (minimum)	30 minutes (minimum)

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY				
Method	Sterilizer cycle	Temperature of exposure	Exposure time	Minimum drying time
Moist heat sterilization	Pre-vacuum	132 <span> </span> °C to 137 <span> </span> °C (270 <span> </span> °F to 277 <span> </span> °F)	4 minutes (minimum)	30 minutes (minimum)

#### H. STORAGE

Sterile instruments should be stored in a clean and dry environment with limited access and under controlled conditions to protect them from dust, moisture, pest and extreme temperatures and humidity. The sterile container should be carefully inspected before use to ensure that it has not been damaged. The instruments are manufactured from non-degradable materials that do not pose any stability problems when stored in the recommended conditions. The responsible user is responsible for ensuring that the instruments, once sterilized, are stored in such a way a to maintain their sterility until the time of use.

#### I. SPECIAL CAUTION: TRANSMISSIBLE SPONGIFORM ENCEPHALOPATHY AGENTS

It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy Agents.

The agents for transmission of Creutzfeldt-Jakob disease are believed to be resistant to normal processes of disinfection and sterilization and therefore the normal processing method of decontamination and sterilization as described above may not be appropriate where Creutzfeldt-Jakob transmission is a risk.

In general, The tissues that come into contact with orthopedic surgical instruments are those of low Transmissible Spongiform Encephalopathy (TSE) infectivity. However, particular precautions should be taken when handling instruments that have been used on known, suspected, or at risk-patients.

#### J. CAUTIONS

- The WasherCap™ Surgical Instrument Set should only be accepted if the original manufacturer’s packaging and labelling are present and arrive intact. If the package has been opened and/or if any component of the Set is damaged, inform your local distributor or Abanza representative for replacement.
- To avoid damaging the instruments that are designed to be turned or screwed must not be hit or blunted. When two devices are intended to be threaded together, ensure that they are fully engaged prior to use.
- Do not use Abanza instruments for any purpose other than their intended use.
- Do not use an instrument that is intended to be used with a specific implant on another implant.
- Some instruments in the set, due to their intended use, have sharp edges. To avoid possible damage to the tendons, patient, and/or users, do not use the instruments for other than the indicated use.
- For the correct use of the countersink instrument, it is required:
  - Disassemble (main body and top element) before carrying out cleaning and disinfection procedures at the reprocessing area, with the help of the countersink tool.
  - Reassemble it for sterilization at the reprocessing area. When assembling, it is important not to mix the parts between the different sizes of countersink (observe the size engraved on each instrument). Check that it has been correctly clamped (when fully threaded) and, if necessary, tighten it with the countersink tool.
  - For further information about the assembly and disassembly of the instruments, please refer to the reprocessing instructions provided by Abanza with the Instrument Set "Instructions for care, cleaning, maintenance and sterilization" (CI-001).

#### K. WARNINGS

- The instruments which compose the WasherCap™ Surgical Instrument Set should be only handled by trauma surgeons who are familiar with the clinical use, surgical technique, proper selection and placement of the device and use of the WasherCap™ associated instrumentation.
- Federal (US) law restricts this device to sale by or on the order of a physician.
- The instruments which compose the WasherCap™ Surgical Instrument Set are supplied non-sterile and must be reprocessed (cleaned, disinfected and sterilized) before use. For this purpose, It is recommended to follow the instructions cited in these Instructions for Use and in the reprocessing instructions provided by Abanza with the Instrument Set "Instructions for care, cleaning, maintenance and sterilization" (CI-001).
- Do not use the instruments for other purposes or surgical procedures other than those indicated in this Instructions for Use.

#### L. DISCLAIMER

Please consult the Instructions For Use (IFU) of the device associated to this surgical instrument set for further information. Please contact your local distributor or Abanza representative for a copy of the IFU of the device associated or the "Instructions for care, cleaning, maintenance and sterilization" (CI-001) of the WasherCap™ surgical Instrument Set.

All surgical instruments which may be damaged, scratched, showing defects, corrosion, cracked joints or discoloration, should not be used and be managed according to hospital protocol rules. Any misuse, modification or substitution of the instruments. Is this recommendation is not followed, Abanza rejects any responsibility derived from any subsequent consequence.

*This document is not subject to control after distribution, and thus, it is the responsibility of the reader to confirm that this is the most current revision or to contact Abanza Tecnomed's customer service to request the latest version.*

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